



MODEL	L IN [mm]	W IN [mm]	HOLE A	B	C IN [mm]	D IN [mm]	SURFACE MATERIAL
BT-1620-SSH	16.25 [412.8]	20.25 [514.4]	1/4-20	1.00 IN	1.63 [41.4]	1.63 [41.4]	LAMINATE WITH STAINLESS STEEL TOP
BT-2024-SSH	20.00 [508]	24.00 [609.6]	1/4-20	1.00 IN	1.50 [38.1]	1.50 [38.1]	
BT-2424-SSH	24.00 [609.6]	24.00 [609.6]	1/4-20	1.00 IN	1.50 [38.1]	1.50 [38.1]	
BT-2430-SSH	24.00 [609.6]	30.00 [762]	1/4-20	1.00 IN	1.50 [38.1]	1.50 [38.1]	
BT-2436-SSH	24.00 [609.6]	36.00 [914.4]	1/4-20	1.00 IN	1.50 [38.1]	1.50 [38.1]	
M-BT-1620-SSH	16.25 [412.8]	20.25 [514.4]	M6	25 mm	1.23 [31.2]	1.27 [32.3]	
M-BT-2024-SSH	20.00 [508]	24.00 [609.6]	M6	25 mm	1.63 [41.4]	1.66 [42.2]	
M-BT-2424-SSH	24.00 [609.6]	24.00 [609.6]	M6	25 mm	1.66 [42.2]	1.66 [42.2]	
M-BT-2430-SSH	24.00 [609.6]	30.00 [762]	M6	25 mm	1.66 [42.2]	1.22 [31.0]	
M-BT-2436-SSH	24.00 [609.6]	36.00 [914.4]	M6	25 mm	1.66 [42.2]	1.27 [32.3]	

- NOTES:
 1. CLEANROOM OPTION HAS THE SAME OUTLINE AND DIMENSIONS.
 2. DIMS ARE IN IN. DIMS IN [] ARE IN mm.